

RECEIVED TEST CARD'S

ScoutPro Quality Control Log – CLIA Waiver

Normal and High Osmolarity Control Solutions are to be tested with each new lot number of Test Cards or when a new shipment of Test Cards is received, even if it is of the same lot number as Test Cards received previously. Record test results below. Follow instructions in the ScoutPro Quick Reference Guide. Call Trukera Medical at (858) 455-6006 if testing falls outside of range.

| Electronic Check Card (ECC) Expected Range: | Normal Control Solution * Optimal Range: | High Control Solution * Optimal Range: |
|--|---|---|
| 331 - 337 mOsm/L | 288 - 305 mOsm/L | 330 - 346 mOsm/L |

| Date | Operator ID | ScoutPro Serial Number | Test Card Lot Number/ Expiration Date | Normal Control Solution | High Control Solution |
|-----------------------|-------------|------------------------|--|--|--|
| Example 10/16/2022 | MSB | 0110060001238 | 00001630 2026/09 | 300 mOsm/L | 341 mOsm/L |
| | | | | <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail | <input checked="" type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |
| | | | | mOsm/L | mOsm/L |
| | | | | <input type="checkbox"/> Pass <input type="checkbox"/> Fail | <input type="checkbox"/> Pass <input type="checkbox"/> Fail |

| *FDA Labeling | Expected | Ranges |
|-------------------------|------------|------------------|
| Normal Control Solution | 297 mOsm/L | 282 – 312 mOsm/L |
| High Control Solution | 338 mOsm/L | 323 – 353 mOsm/L |